

Product Name :
Emergency Safety Shower and Eye Wash

Product Code :
LIM-CAT-L0021-00018



Description :

Emergency Safety Shower and Eye Wash

Technical Specification :

Emergency Safety Shower and Eye Wash
Emergency Self-Draining Safety Shower and Eye Wash which meet the following general specifications:
Weight: 16kgs
Dimension (W*D*H): 580mmx750mmx 2300mm
Material of construction:
Pipework and fittings SUS 316L and fittings
Operating valves SUS316L, two piece stay open valves
Shower Spray Outlet Machine molded nylon 6 shower nozzle
Eye/Face Wash Bowl/Lid SUS316L stainless steel
Eye Wash Diffusers Chrome plated brass with SUS316L mesh
Standard Finish: Polished brush finish
Service:
Water Inlet: 1 1/4" Female BSP
Design Operating Pressure: 2.0 to 6 BAR
Test Pressure: up to 10.3 BAR
Flow at Operating Pressure Shower: 76liters/minute
Eye/face wash: 12 liters/minute
Environmental Operating Temperature:
Ambient operating temperature: 5°C

Ambient Operating Temperature: 45°C
Operational: shower spray operation: Hand Lever.

Emergency Safety Shower and Eye Wash, Emergency Safety Shower and Eye Wash Bulk Suppliers, Emergency Safety Shower and Eye Wash Tools, Emergency Safety Shower and Eye Wash Emergency Safety Shower and Eye Washes, Emergency Safety Shower and Eye Wash Manufacturers, Emergency Safety Shower and Eye Wash Suppliers from India, China, Kenya.



Laboratory instruments manufacturers India