

Product Name :
Strips Cups

Product Code :
LIM-CAT-L0042-00024



Description :

Strips Cups

Technical Specification :

Strips cups
Plastic green (standard in the dairy industry)
The strip cup is suitable for checking the first jets of milk on flakes, for the detection of mastitis.
Including strainer
Result immediately visible
Mastitis test
Anti-splash screen
Sieve can be easily unclipped and cleaned
Robust plastic.

Strips Cups, Strips Cups Bulk Suppliers, Strips Cups Tools, Strips Cups Strips Cupss, Strips Cups Manufacturers, Strips Cups Suppliers from India, China, Kenya.

