

Product Name :
Ear Notches

Product Code :
LIM-CAT-L0042-00029



Description :

Ear Notches

Technical Specification :

Ear notches

Standard model made of metal.

For larger pigs, a notcher that makes notches 1/2 inch deep is recommended

Ear-notching helps identify a pig's litter and which one of the litter it is, giving each pig a unique identity number.

Notches are placed in one of five locations in the pig's right ear to show the litter number and in one of three locations in the left ear to show the individual pig number.

Ear Notches, Ear Notches Bulk Suppliers, Ear Notches Tools, Ear Notches Ear Notchess, Ear Notches Manufacturers, Ear Notches Suppliers from India, China, Kenya.