

**Product Name :**  
Dehorning Saw

**Product Code :**  
LIM-CAT-L0042-00044



**Description :**

Dehorning Saw

**Technical Specification :**

Dehorning Saw

Length 36 cm with wooden handle.

In this method of dehorning, a hand saw is used to remove the horn along with a ½-inch ring of skin around the horn base.

Advantages: Can be used in older calves with large horns.

Disadvantages: Risk of infection and death from blood loss; may lead to set-backs in overall health and performance, including weight loss; animal is monitored for several days afterward.

Dehorning Saw, Dehorning Saw Bulk Suppliers, Dehorning Saw Tools, Dehorning Saw Dehorning Saws, Dehorning Saw Manufacturers, Dehorning Saw Suppliers from India, China, Kenya.

