

**Product Name :**  
Dehorning Wire

**Product Code :**  
LIM-CAT-L0042-00045



**Description :**

Dehorning Wire

**Technical Specification :**

Dehorning Wire

Wire type consisting of 2 handles and 40ft/12 mtrs special cutting wire.

The Dehorning Wire is a length of stainless steel saw wire used to remove mature horns from adult cattle, rams and goats.

The wire cuts through the horns quickly and without any slippage.

Use with the handles supplied.

Dehorning Wire, Dehorning Wire Bulk Suppliers, Dehorning Wire Tools, Dehorning Wire Dehorning Wires, Dehorning Wire Manufacturers, Dehorning Wire Suppliers from India, China, Kenya.