

**Product Name :**  
Soxhlet Extractor

**Product Code :**  
LIM-CAT-L0046-00012



**Description :**

Soxhlet Extractor

**Technical Specification :**

Soxhlets are made of low extractable borosilicate glass and come complete with extractor, condenser, and flask. An extraction thimble is required to support the sample within the extractor. Place a soil sample in the thimble and heat the solvent contained in the flask. As the solvent vapors rise (passing through the extraction thimble), they enter the water-cooled condenser, and reliquify. When the liquid level in the extractor reaches the top of the siphon tube, siphoning action returns the extract-enriched solvent to the flask.