

**Product Name :**  
Gas/Liquid Relative Permeameter

**Product Code :**  
LIM-CAT-L0046-00014



**Description :**

Gas/Liquid Relative Permeameter

**Technical Specification :**

Pressure Gauge – 0.00 to 1.00 ATM

Golden Core Holder to hold different sizes of core ½” to 1 ½” through rubber ring.

Digital Temperature indicator attached with core holder

A Flow meter (Rota meter) of capacity 0 to 250 cc/sec with regulator holding arrangement.

Regulator to operate different pressure through a gas pressure line attached with pressure gauge.