

**Product Name :**  
Variable Speed Electronic Rotational Viscometer

**Product Code :**  
LIM-CAT-L0047-00001



**Description :**

Variable Speed Electronic Rotational Viscometer

**Technical Specification :**

**Anti-Slip Design** - The handle is made of a matte process, and the frame is made of high-strength magnesium alloy and 304 stainless steel to achieve anti-slip effect, prevent the error caused by hand slip. The fuselage uses an all-metal structure to extend the life of the machine.

**LCD Display** - The viscometer is equipped with a high-definition LCD display with 6 operating buttons, which can be operated according to the actual use and get clear data display.

**4 Speed Adjustment** - The viscometer offers an adjustable speed, 6 rpm, 12 rpm, 30 rpm, 60rpm revolution/min, which can be adjusted according to the specific needs.

**Wide Application** - The digital viscometer is widely used in fluid, semi-fluid production industries and scientific research units such as oils, paints, coatings, slurries, textiles, foods, drugs, pharmaceuticals, grease, adhesives and cosmetics etc.

**Features:**

Liquid crystal display with six button control.

Speed adjustable to meet your various requirements.

Grinding process to achieve a non-slip effect.

A set of supporting tools can provide a better user experience.



## Laboratory instruments manufacturers India