

Product Name :
Lab Electric Stirrers For Mixing Solutions

Product Code :
LIM-CAT-L0047-00012



Description :

Lab Electric Stirrers For Mixing Solutions

Technical Specification :

Lab electric stirrers for mixing solutions
20 liters capacity
Advanced Stirrers with microprocessor technology
Speed range 50 – 2200 rpm
Monitoring of set & actual speed
Constant speed even with change of load / Voltage
High torque even at low speeds
Accurate stirring speed control
Maintenance free Brushless DC motor ensures high torque
Easy-to-read backlit LCD for display of set & actual parameters
Display of torque / speed
Last parameter recall, ideal for repetitive processes
Motor protection in case of continuous overload
Accepts shafts up to 13mm Ø
Stopper for protection of vessel & shaft
Special base design for topple free operation.



Laboratory instruments manufacturers India