

Product Name :
Mud Balance

Product Code :
LIM-CAT-L0047-00014



Description :

Mud Balance

Technical Specification :

Mud Balance is Designed to find out specific gravities of semi liquids like mud and other liquids.
It has a stainless steel beam calibrated.
A stainless steel cup with lid and overflow vent is fitted on one side of the beam.
A counter weights with cursor slides over the graduated scale.
The beam has a knife-edge at centre which rests in a fulcrum fitted in the stand.
Leveling screws and spirit level are fitted to the stand.
Supplied complete with box.