


Product Name :
Conductivity Meter

Product Code :
LIM-CAT-L0048-00013



Description :

Conductivity Meter

Technical Specification :

Digital Conductivity Meters are reliable and accurate test instruments for measurement of Conductivity of aqueous solutions. They measure Conductivity in five ranges.

The Resolution is 0.1S/Cm in the lowest range. The cell constant is flashed on a digital display & can be adjusted from the front panel. The instruments have the Check facility to calibrate the instrument.

These instruments are ideal for monitoring salt contents in water and soluble salts in soils. They are extremely useful instruments for agriculture and soil analysis labs, swimming pools, fertilizer plants, petroleum refineries, textile plants etc.

Features:

Highly Stable and Accurate

Auto Temperature Compensation Facility

3 Digit Display.
