



Email :
sales@laboratoryinstrumentsmanufacturers.com

Product Name :
Sur Charge Weight 2.5kg

Product Code :
LIM-CAT-L0049-00021



Description :

Sur Charge Weight 2.5kg

Technical Specification :

Sur Charge Weight 2.5kg
Surcharge Weights are used with CBR or LBR testing in the application of surcharge forces.
They are made of plated steel and available in annular or slotted models, with a 5.9 in (149mm) diameter.

Features:

Plated steel construction
Available in annular or slotted design



Laboratory instruments manufacturers India