

**Product Name :**  
Electronic Weigh Balance 600gm

**Product Code :**  
LIM-CAT-L0049-00025



**Description :**

Electronic Weigh Balance 600gm

**Technical Specification :**

Electronic weigh balance 600gm  
Body high grade engineered plastic with satin finish is easy to clean with the help of damp cloth. 110 mm round Pan for easy loading of Tare container or item directly on the pan is made of stainless steel.  
Runs on AAA size alkaline battery pair as well as 220 VAC through adapter (Both options provided with the scale) It saves on running cost by avoiding replacement of batteries.  
Simple single touch calibration is easy to perform. Multiple weighing modes available to convert units from g to oz or to lb during weighing.