

**Product Name :**  
Electronic Weigh Balance 5kg

**Product Code :**  
LIM-CAT-L0049-00026



**Description :**

Electronic Weigh Balance 5kg

**Technical Specification :**

Electronic Weigh Balance 5kg  
Large font LCD display with backlight, clear reading.  
Stainless Steel Platter  
Gram, LB, OZ, CT. Unit  
Simple 5 keys panel for easy operation  
Adjustable feet  
Overload Protection  
Below Balance Weighing Facility  
Checking Weighing, Piece counting function  
External power supply adapts and battery inside.

