

Product Name :
Sieve Shaker Machine

Product Code :
LIM-CAT-L0060-00002



Description :

Sieve Shaker Machine

Technical Specification :

The analytical sieve shakers of the are used in research & development, quality control of raw materials, interim and finished products as well as in production monitoring.

The controllable electromagnetic drive offers an optimal adaption for every product.

Sharp fractions are obtained even after very short sieving times.

Mechanism: Consists of Sieving mechanism, to accommodate 10 nos. of sieves of 200mm Dia. & 35mm Ht

Supply specifications to adapt equipment accordingly (Default: 230v AC, 50/60Hz).

Set of sieves 200 mm. Dia& 35 mm. Height of each standard sieve.

Sieve set consists of sieves with ISS No. in microns 1700, 850, 600, 425, 300, 212, 150, 125, 106, 75, 53 with lid & receiving pan

Width/ depth/ height, mm 470*330*630

Weight 44kg.

Range 100gsm.

