

**Product Name :**  
Defibrillator Kit

**Product Code :**  
LIM-CAT-L0072-00004



**Description :**

Defibrillator Kit

**Technical Specification :**

Defibrillator Kit

Fully automatic defibrillator - Fully automatic delivery of the shock on demand, without user intervention.

Self-tests of the battery, internal circuits and the presence and function of the defibrillator electrodes.

Configurable voice prompts.

Customizable voice support in 3 levels between very detailed instructions for lay users and short sentences for experienced rescuers.

Speech and text prompts tailored to the pace of the user.

Backlit LCD display.

Circumferential rubber coating IP55 degree of protection against the ingress of dust and water.