

Product Name :
Tripod accessories for Performing Experiments in Physics

Product Code :
LIM-CAT-L0018-00002



Description :

Tripod accessories for Performing Experiments in Physics

Technical Specification :

Tripod accessories for performing experiments in physics

Tripod stand with base:

Height: 750mm

Diameter: 12mm

Base made of cast iron

Possibility of removing of tripod rod

Dimensions of the base: (210*130) mm

Tripod rod:

Material: stainless steel

M10 final processing

Length: 1000mm

Diameter: 12mm

Tripod accessories for Performing Experiments in Physics, Tripod accessories for Performing Experiments in Physics Bulk Suppliers, Tripod accessories for Performing Experiments in Physics Tools, Tripod accessories for Performing Experiments in Physics Tripod accessories for Performing Experiments in Physics, Tripod accessories for Performing Experiments in Physics Manufacturers, Tripod accessories for Performing Experiments in Physics Suppliers from India, China, Kenya.



Laboratory instruments manufacturers India