

**Product Name :**  
Laboratory Microwave Oven

**Product Code :**  
LIM-CAT-L0082-00008



**Description :**

Laboratory Microwave Oven

**Technical Specification :**

Oven is fitted with heating elements on all three sides, ensuring uniform temperature distribution throughout the working chamber  
Sealed and enclosed heaters to working chamber for minimum heat loss  
Air heaters placed in the air path created at the back for uniform heating and uniform heat distribution Temp. Range.  
Special heavy-duty Stainless-steel lock and door hinges with spring and roller mechanism  
Stainless steel crimped wire mesh Shelves to withstand heavy load and to offer minimum resistance to air circulation for better heat distribution  
Provision to add more shelves, as per requirement of test - Adjustable Shelves to place test samples or Glass wares of varying dimensions  
Adjustable air ventilation mounted at top  
Operation on 220/230 Volt, 50 Hz, Single phase, AC Supply  
With built-in forced air circulation arrangement for uniform heat distribution inside the working chamber  
Complete with pilot indicating lamps, 2 meter cord wire and power plug.

---

## Laboratory instruments manufacturers India

