

**Product Name :**  
RF Signer Generator Trainer

**Product Code :**  
LIM-CAT-L0084-00010



**Description :**

RF Signer Generator Trainer

**Technical Specification :**

Wide frequency range: 10-600MHz.

Internal/external modulation (AM/FM) Fine control over Frequency / Amplitude range Bright 4 1/2 digit LED/16X2 LCD display for freq.

Two models- manually operated or micro-controller based.

Useful for testing and measurement in RF laboratory.

Optionally PC GUI based control using USB IO Module (UCOM).