

**Product Name :**  
Set for performing experiments in the field of Light box and Optics

**Product Code :**  
LIM-CAT-L0018-00010



**Description :**

Set for performing experiments in the field of Light box and Optics

**Technical Specification :**

Set for performing experiments in the field of Light box and Optics.

Set of 5 acrylic prisms, and each has rectangle, semicircle, Equilateral triangle  $60^\circ \times 60^\circ \times 60^\circ$ , triangle  $90^\circ \times 60^\circ \times 30^\circ$  i Right angle triangle  $90^\circ \times 45^\circ \times 45^\circ$ . All prisms are without defects and all have sharp edges.

Set of 3 acrylic lenses: one is double concave, one is double convex (both are having the same radius of curvature) and one is thick double convex lens.

Set of 3 mirrors: polished, reflexing surfaces and freestanding. Every set has a parabolic, semicircular and flat glass.

Set of 2 separated plates, black: one is with 1 i 3 narrow notch, and second is with 4 narrow notch i 1 wide notch.

Set of 8 filters in colours: 3 primary colours (red, blue i green) i 5 sekundary colours (cyan, violet, orange, yellow i magenta).

Set of 8-colored plates: red, blue, green, violet, orange, cyan, yellow i pink.

Set delivers with a backup bulb in the styrofoam box, for longer time of storing.

Set for performing experiments in the field of Light box and Optics, Set for performing experiments in the field of Light box and Optics Bulk Suppliers, Set for performing experiments in the field of Light box and Optics Tools, Set for performing experiments in the field of Light box and Optics Set for performing experiments in the field of Light box and Opticss, Set for performing experiments in the field of Light box and Optics Manufacturers, Set for performing experiments in the field of Light box and Optics Suppliers from India, China, Kenya.



## Laboratory instruments manufacturers India