

**Product Name :**  
Optical Bench

**Product Code :**  
LIM-CAT-L0018-00011



**Description :**

Optical Bench

**Technical Specification :**

Optical bench

Set made of Aluminium, Includes a leveling screw,

Scale is divided in millimeters,

It has antitilt guides, with an index mark for mounting various accessories,

Height of accessories can be adjusted.

Optical bench delivers with the following accessories:

Halogen lamp,

Flexible cables,

Round flat glass,

Lens holder

Round black screen, diaphragm, ground glass screen, diffraction grating holder, large slits, Jolly's photometer, pin support and receiving screen.

Optical Bench, Optical Bench Bulk Suppliers, Optical Bench Tools, Optical Bench Optical Benches, Optical Bench Manufacturers, Optical Bench Suppliers from India, China, Kenya.

---

## Laboratory instruments manufacturers India

